



Rev 1: September 2018

FSN Ref: DNB.444.489.2023

FSCA Ref: DNB.444.489.2023

VAIM Co., Ltd.

H O: 103, Okcheon Medical Device Incubation Center, 79, Uryodanji-gil,  
Okcheon-eup, Okcheon-gun, Chungcheongbuk-do, 29055, K  
T.+82.43.731.6456 F.+82.43.731.6455

Daejeon: 14, Techno 3-ro, Yuseong-gu, Daejeon, 34012, Korea  
T.+82.42.935.6456 F.+82.42.935.6455

Date: 02.06.2024

## **Urgent Field Safety Notice**

### **Juvelook**

**For Attention of\*:** Identify either by name or role who needs to be aware of the hazard and/or take action. If this is multiple recipients then include full list.

<b>Contact details of local representative (name, e-mail, telephone, address etc.)*</b>
Beautyeuropa.eu Sp.z o.o.
Filip Skrzypiec – President Włodzimierz Lewitin – National Sales Manager
Skowronków 10 I 58-500 Jelenia Góra, Poland
<a href="mailto:w.lewitin@beautyeuropa.eu">w.lewitin@beautyeuropa.eu</a>
+48 518 225 625

**Urgent Field Safety Notice (FSN)**

**Juvelook**

**Risk addressed by FSN**

<b>1. Information on Affected Devices*</b>	
1	<b>1. Device Type(s)*</b>
.	Soft Tissue injectable Implant / Class III
1	<b>2. Commercial name(s)</b>
.	Juvelook
1	<b>3. Unique Device Identifier(s) (UDI-DI)</b>
.	(01)08800051300102
1	<b>4. Primary clinical purpose of device(s)*</b>
.	Device is intended to be used for intradermal and subcutaneous implantation to fill skindepressions as to restore and enhance soft tissue volume.
1	<b>5. Device Model/Catalogue/part number(s)*</b>
.	Juvelook
1	<b>6. Software version</b>
.	N/A
1	<b>7. Affected serial or lot number range</b>
.	VJ220511
1	<b>8. Associated devices</b>
.	N/A


<b>2 Reason for Field Safety Corrective Action (FSCA)*</b>	
2	<b>1. Description of the product problem*</b>
.	The adverse reaction occurred after administration of Juvelook (Lot no. VJ220511) on 3rd July 2023. Juvelook was applied to the eye area. After administration of the preparation, large thickenings (deposits, lumps) occurred.  According to Polish CA opinion, in Juvelook IFU there is no clear description of area where medical device should be injected or clear delimitation of safe and non-safe areas. And, incident was caused by lack of clear relevant information on device's IFU about restrictions for application of Juvelook into eyelids.
2	<b>2. Hazard giving rise to the FSCA*</b>
.	The product must be handled and used by a plastic surgeon, dermatologist or a trained doctor with special knowledge and skills after fully understanding the product itself and its directions for use.
2	<b>3. Probability of problem arising</b>
.	None.
2	<b>4. Predicted risk to patient/users</b>
.	None.
2	<b>5. Further information to help characterise the problem</b>
.	Juvelook should not be used in periorbital area (eyebrow, eye spacing, eyelashes, upper eyelid, and lacrimal punctum)

	<p><b>【Warnings】</b>          Serious side effects such as loss of vision may occur when injected into blood vessels. Therefore, it should not be used in periorbital area (eyebrow, eye spacing, eyelashes, upper eyelid, lacrimal punctum).</p> <p>Do not inject products into shallow skin layers in periorbital area (eyebrow, eye spacing, eyelashes, upper eyelid, lacrimal punctum) to prevent papules or nodules. Papules or nodules can be caused by inappropriate procedures (injection into skin surface, an excessive amount of injection). Massaging the injection site to distribute the product evenly can minimize the incidence of papules and nodules.</p> <p><b>【General Precautions】</b>          Care should be taken not to inject thinner skin areas, such as periorbital area, because of the high risk of papules or nodules.</p>
2	<p><b>6. Background on Issue</b></p> <p>We've only contracted with Beautyeurope.eu Sp.z o.o., our sole official representative and exclusive distributor, of Lensina and Juvelook products in Poland. The patient purchased the Juvelook through Lipstore which is not our distribution channel in Poland and was not injected by well-trained surgeon, dermatologist or a trained doctor for our product.</p>
2	<p><b>7. Other information relevant to FSCA</b></p> <p>Product purchases must be made through established distribution channels.</p>

<b>3. Type of Action to mitigate the risk*</b>			
<b>3.</b>	<p><b>1. Action To Be Taken by the User*</b></p> <p> <input type="checkbox"/> Identify Device    <input type="checkbox"/> Quarantine Device    <input type="checkbox"/> Return Device    <input type="checkbox"/> Destroy Device  <input type="checkbox"/> On-site device modification/inspection  <input type="checkbox"/> Follow patient management recommendations  <input checked="" type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU)  <input type="checkbox"/> Other                      <input type="checkbox"/> None         </p>		
<b>3.</b>	<table border="1" style="width: 100%;"> <tr> <td style="width: 30%;"><b>2. By when should the action be completed?</b></td> <td style="width: 70%;">Specify where critical to patient/end user safety 1Q, 2025</td> </tr> </table>	<b>2. By when should the action be completed?</b>	Specify where critical to patient/end user safety 1Q, 2025
<b>2. By when should the action be completed?</b>	Specify where critical to patient/end user safety 1Q, 2025		
<b>3.</b>	<p><b>3. Particular considerations for:</b> Choose an item.</p> <p>Is follow-up of patients or review of patients' previous results recommended? Choose an item.</p> <p>Provide further details of patient-level follow-up if required or a justification why none is required</p>		
<b>3.</b>	<table border="1" style="width: 100%;"> <tr> <td style="width: 70%;"><b>4. Is customer Reply Required? *</b> (If yes, form attached specifying deadline for return)</td> <td style="width: 30%; text-align: center;">No</td> </tr> </table>	<b>4. Is customer Reply Required? *</b> (If yes, form attached specifying deadline for return)	No
<b>4. Is customer Reply Required? *</b> (If yes, form attached specifying deadline for return)	No		

3.	<p><b>5. Action Being Taken by the Manufacturer</b></p> <p> <input type="checkbox"/> Product Removal                      <input type="checkbox"/> On-site device modification/inspection  <input type="checkbox"/> Software upgrade                      <input checked="" type="checkbox"/> IFU or labelling change  <input type="checkbox"/> Other    <input type="checkbox"/> None </p> <p>Provide further details of the action(s) identified.  In the amended IFU, there will be clear direction that periorbital area includes eyebrow, eye spacing, eyelashes, upper eyelid, and lacrimal punctum and Juvelook shouldn't be injected in periorbital area in <b>【Warnings】</b>and<b>【General Precautions】</b>of Juvelook IFU.</p>	
3	6. By when should the action be completed?	<p>1Q, 2025</p> <p>*Before the implementation of new Juvelook IFU, we will provide users with safety notes stating corrected warnings and general precautions.</p>
3.	7. Is the FSN required to be communicated to the patient /lay user?	No
3	<p>8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?</p> <p>Choose an item.                      Choose an item.</p>	

<b>4. General Information*</b>	
4.	1. FSN Type* <span style="float: right;">New</span>
4.	2. For updated FSN, reference number and date of previous FSN <span style="float: right;">Provide reference and date of previous FSN if relevant</span>
4.	3. For Updated FSN, key new information as follows: Summarise any key difference in devices affected and/or action to be taken.
4.	4. Further advice or information already expected in follow-up FSN? * <span style="float: right;">No</span>
4	5. If follow-up FSN expected, what is the further advice expected to relate to: Eg patient management, device modifications etc
4	6. Anticipated timescale for follow-up FSN <span style="float: right;">For provision of updated advice.</span>
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)
	a. Company Name <span style="float: right;">Only necessary if not evident on letter-head.</span>
	b. Address <span style="float: right;">Only necessary if not evident on letter-head.</span>
	c. Website address <span style="float: right;">Only necessary if not evident on letter-head.</span>
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. * Yes
4.	9. List of attachments/appendices: <span style="float: right;">N/A</span>
4.	10. Name/Signature <span style="float: right;">Seungwook Lee Head of Regulatory Affairs</span>



<b>Transmission of this Field Safety Notice</b>	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback..*</p>

Note: Fields indicated by \* are considered necessary for all FSNs. Others are optional.